

# DISPOSABLE PROTECTIVE FACE MASKS

## FOCUS PRODUCT

### REFERENCE

**Blue mask:**

1677742680

**Black mask:**

1677742780



→ Individual disposable masks for repairers.

→ Good Quality worth cost.

→ For enhanced hygiene: tamper-evident label and protective bag on the box.

### Precautions for use

Before putting the mask, wash your hands with water and soap or with an hydro alcoholic solution.

- When wearing a mask, do not touch it or displace it.
- Once you touch a used mask, wash your hands with water and soap or with an hydro alcoholic solution.
- If you need to drink or eat something, change your mask.
- When a mask is moistened, replace it with a new one, do not re-use, single-use masks.
- To remove the mask :
  - Remove it from the ears, do not touch the front of the mask,
  - Throw the mask in a closed trash once,
  - Wash your hands with water and soap or with an hydro alcoholic solution,
- Do not use if you are allergic to non-woven materials.
- Do not use if the packaging is damaged or open.
- Disposable mask: must be destroyed after use, no washing and no re-use.
- Keep in a clean and ventilated area.
- Not intended for use by healthcare professionals in an operating room or in other medical settings with similar requirements.

### TECHNICAL FEATURE

- **Type II compliant:**
  - Bacterial Filtration Efficiency (BFE) > 99,98 %,
  - Differential Pressure : 22,57 Pa/cm<sup>2</sup>,
  - Microbal load : 3,69 CFU/g,
- 50 masks per box

### REGULATION

- Conformity declaration and test report are attached to the product.
- CE certification
- EN 14683 type 1 standard

### GOOD TO KNOW

- With a secure bag in the box, the **masks can be offered for sale** to garage customers.
- **The box is placed on the counter**, the receptionist distributes a mask with gloves.

